NiTiShaper

Dispositif médical de classe IIa

F.F.D.M. PNEUMAT

78-80 avenue de la Prospective

18020 BOURGES Cedex

CHER / CENTRE / France

Tél. : 02 48 23 27 37

thomas.dental@ffdm-pneumat.com

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INDICATIONS

For dentin removal and root canal shaping.

CONTRE-INDICATIONS

This product contains Nickel: do not use NiTiShaper on patients who are sensitive or allergic to Nickel. Do not use NiTiShaper if an instrument is damaged or contaminated.

PREAMBULE

NiTiShaper is aimed to be used by professionals of dental art only during root canal treatment. Professionals learn how to use these instruments and the related medical procedures during their training.

It is recommended to use NiTiShaper on an extracted tooth or endo training block prior to using it on a patient.

DESCRIPTION

NiTiShaper is a kit of 3 endodontic instruments. The blades of the instruments are made of a nickel-titanium alloy. The part ensuring the connection with the contra angle is manufactured according to the ISO standard 1797 (diameter of the contra angle shank: $\varnothing 2.35$ / -0.016 mm). Supplied with a silicone stopper on the blade (silicone USP class 6).

PACKAGING

Quantity per pack: 3 units or 6 units, identical or different (assortment). Blister non sterile.

LIFE SPAN

The life span is estimated at 3 years from the date of manufacturing but can vary depending on the use (cutting, conditions), on the frequency of use, and on the quantity of sterilization.

STOCKAGE

Store NiTiShaper in a clean, dry and dust free area. Apply the "first-in first-out" method for the stock management.

PROTOCOL OF USE

Before starting the clinical procedure, X-ray examination must be made to evaluate the root canal anatomy (width, length, curvature).

Coronal access

Before using NiTiShaper, create an access cavity to expose all the canal orifices, to allow easy viewing, and straight line instrumentation and to ensure a good irrigation during all the treatment. For this operation, we recommend to use Gates, RA steel Round burs, Endopener from THOMAS.

Glide path

Create a glide pathway by using a manual canal probe ISO 12 or a K File ISO 10. Irrigate and introduce the K file with a reciprocating movement, into the root canal.

This step is very important to establish the working length and patency of root canals with the small size K File in stainless steel. K-Files have excellent tactile sensation. The stiffness of hand steel K Files helps to negotiate blockages and calcifications.

This step will facilitate the penetration of the other instruments.

This step must be realized between each NiTiShaper (T1, T2, T3).

STEP 1 : T1 : SECURE OF THE ROOT CANAL

Taper : 4%
Tip Ø : 0.18mm
Length of the blade : 25mm
Color code : white
Depth marks : 20, 22mm
Total length : 37mm

Introduce in rotation, the first NiTiShaper instrument "T1" – white strip. Apply a slight pressure with vertical filling movements. You can reach the working length.

- Rotation speed: 300 to 500 rpm
- Torque: 3N.cm to 4N.cm

Clean and irrigate the root canal (see Glide path above).

STEP 2 : T2 : SHAPING

Taper : 6%
Tip Ø : 0.20mm
Length of the blade : 25mm
Color code : red
Depth marks : 18, 19, 20, 22mm
Total length : 37mm

Introduce in rotation, the second NiTiShaper instrument "T2" – red strip. Apply a slight pressure with vertical filling movements. You can reach the working length.

For curved root canal, we recommend to introduce first "T2" at 2/3 of the working length, then the K File.

Then, you can reach the working length with T2".

- Rotation speed: 300 to 500 rpm
- Torque: 3N.cm to 4N.cm

Clean and irrigate the root canal (see Glide path above).

ETAPE 3 : T3 : FINITION

Taper : 6%
Tip Ø : 0.23mm
Length of the blade : 25mm
Color code : black
Depth marks : 18, 19, 20, 22mm
Total length : 37mm

Introduce in rotation, the third NiTiShaper instrument "T3" – black strip. Apply a slight pressure with vertical filling movements. You can reach the working length.

- Rotation speed: 300 to 500 rpm
- Torque: 3N.cm to 4N.cm

Clean and irrigate the root canal (see Glide path above).

MAINTENANCE and STERILIZATION

- It is essential to follow the steps of cleaning, disinfection and sterilization before the first use or for re-use. NiTiShaper is supplied nonsterile.
- Open the alveolus and make the sterilization in a thermo welded bag (check that bags comply with EN ISO 11607-1 standard). Autoclave 134°C during 18 minutes.
- It is the user's responsibility to sterilize products before first use and re-use.
- Dispose of the product if damaged or decreased in efficiency. It is the user's responsibility to evaluate the level of wear after the first use. In order to avoid any risk, we recommend a single use.
- For your safety, we recommend to use Individual Protective Equipment (IPE).
- Use NiTiShaper carefully to avoid puncturing finger.
- Only use disinfectants with proven effectiveness (EU labeling, FDA approved...).
- Follow strictly the instruction manual of the manufacturer of the medical cleaning product.
- When several instruments are sterilized in a single autoclave cycle, ensure that the maximum load of the autoclave is not exceeded.

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LOT LOT - Ce symbole est accompagné du code du lot du fabricant / The symbol is accompanied by the manufacturer's lot code.

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