

# Product Catalogue




# incidental




MADE IN  
TÜRKİYE

**ENDOART®**  
FINAL TOUCH



**ruby | dent**





# INDEX

INTRODUCTION	5	RubyComp P Posterior Composite	25
ENDODONTICS		RubyFlow Flowable Composite	25
EndoArt Action Rotary Files	7	RubyBond Adhesive System	26
EndoArt Smart Rotary Files	8	RubySEBond Self Etching Adhesive System	26
EndoArt Touch Rotary Files	9	RubyEtch Etching Gel	27
EndoArt Expert Reciprocating Files	11	RubyEtch HF Porcelain Etch	27
EndoArt Pedo Rotary Files	12	RubyLiner with Glass Ionomer	28
EndoArt Retreatment Files	13	RubyCal LC Calcium Hydroxide	28
EndoArt Path Rotary File	13	Ruby Pit and Fissure Sealant	29
EndoArt K Files	14	FlexDam Gingiva Barrier	29
EndoArt H Files	14	Rubydent Curing Light	30
EndoArt Spreader	14	RubyShine Prophylaxis Paste	31
EndoArt C Files	14	RubyCap Capsule Amalgam	32
EndoArt Endomotor with Apex Locator	15	RubyCap NG Capsule Amalgam	32
EndoArt Apex Locator	15	RubyCap High Capsule Amalgam	33
EndoArt Gutta Percha Points	16	RubyGam Powder Amalgam	34
EndoArt Absorbent Paper Points	16	RubyGam High Powder Amalgam	34
EndoArt Fiber Post	17		
RubyCal Calcium Hydroxide Paste With Barium Sulphate	18	POLISHING & FINISHING	
RubyCal Calcium Hydroxide Paste With Iodoform	18	RubyPlaton Polishing Discs	36
RubyCal Calcium Hydroxide Powder & Liquid	18	RubyPlaton Polishing Strips	36
RubyCal Calcium Hydroxide Powder	18		
RubyZoe Zinc Oxide Eugenol Cement	19	WEDGES	
RubyZoe Eugenol	19	RubyPlaton Wooden Wedges	37
RubyZoe Zinc Oxide Powder	19	RubyPlaton Transparent Wedges	37
RubyEdta EDTA Gel	20		
EndoArt EDTA Solution	20	ARTICULATING PAPERS	
EndoArt Chlorex Gel	20	RubyMark Articulating Paper Straight	38
EndoArt Chlorex	20	RubyMark Articulating Paper Horseshoe	39
EndoArt Sodium Hypochloride	21	RubyMark Articulating Paper 200µ Straight	39
EndoArt Citric Acid	21		
Ruby Iodoform Powder	21	AUXILIARY	
Ruby Gutta Solver	21	RubyCrownRemover	40
		RubyPlaton Cotton Roll	40
RESTORATIVES			
RubyFill Nano Composite	24	STUDY MODELS	
RubyComp Nano Composite	24	Arma Dental Study Model	41
		Arma Dental Phantom Teeth	41



# About us

## **Company**

Inci Dental was founded 1994 in Istanbul. Since then the company produce a variety of consumable materials. Since 1994, the concept of quality is established, always giving priority to customer satisfaction, the production process by making significant investments in R & D activities, cooperate with universities and Scientific and Technological Research Council of Turkey, made production as newest and unique products, in the sector so has achieved a leading place.

The company is only producer of variety materials in Turkey. Inci Dental is manufacturer of Dental amalgam filling material, dental composites filling material, endodontic rotary files, fiber posts, polishing discs, articulating paper and etc. We export our products more than 60 countries in the World.

Our quality assurance approve by European Union with the certificates of ISO 9001, ISO 13485 and the EC.

## **Vision**

To be a company which is the first choice of domestic and foreign business partners.

In this context, the sale of the products that we have done all, are the first brands that come to mind of the dentists.

## **Mission**

Inci dental, in dental sector;

Who makes the difference

Who find, research and use the right technologies in all over the world

Who adapt quickly to technological developments

Who produce best quality

Who lead the innovation

A company which holds the highest level of employee, domestic and foreign business partners satisfaction.



**ENDOART**<sup>®</sup>  
FINAL TOUCH



**Heat Treated Control Memory  
Rotary Files**



**ACTION**

SX	S1	S2	F1	F2	F3	F4	F5		SX	S1	S2	F1	F2	F3	F4	F5	
Starter		Shaping		Finishing					Starter		Shaping		Finishing				
Torque: 3.0 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm		Torque: 2.6 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 1.8 Ncm 300 Rpm	Torque: 2.6 Ncm 300 Rpm

↑ Excellent Cutting Efficiency     ↻ Maximum File Flexibility

EndoArt Action rotary files allow you to move comfortably in the canal thanks to the special structure of its surface, while its variable taper is designed to evacuate debris from the canal. The developed special heat treatment technology ensures that the Ni-Ti wire has a much more flexible and durable structure and also turns it into a Gold or Blue color. EndoArt Action Gold files have a higher sharpness efficiency and higher fracture resistance than standard Ni-Ti files. EndoArt Action Blue files allows pre-bending before entering the canal thanks to its high flexibility efficiency with special control memory and it has a much higher fracture resistance than standard Ni-Ti files. It adapts to the shape of the canal in curvature structures.

Action Kit    SX - S1 - S2 - F1 - F2 - F3

**Action Gold Rotary Files**

- Rotation Files
- Variable Taper
- Better Cutting Efficiency
- x2 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box



**Action Blue Rotary Files**

- Rotation Files
- Variable Taper
- Better Flexibility
- X4 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box

**Wire** : Gold  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files

**Wire** : Blue  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files



**WhatsApp**  
Get information!

SMART



10/04	15/04	20/04	25/04	30/04	35/04	40/04	50/04	10/04	15/04	20/04	25/04	30/04	35/04	40/04	50/04																
<b>04 Taper</b>		Starter		Shaping / Finishing																											
Torque: 3.0 Ncm 300 Rpm		Torque: 2.0 Ncm 300 Rpm		Torque: 2.0 Ncm 300 Rpm		Torque: 2.0 Ncm 300 Rpm		Torque: 2.0 Ncm 300 Rpm		Torque: 2.0 Ncm 300 Rpm		Torque: 2.0 Ncm 300 Rpm		Torque: 2.0 Ncm 300 Rpm		<b>04 Taper</b>		Starter		Shaping / Finishing											
Torque: 3.0 Ncm 300 Rpm		Torque: 1.6 Ncm 300 Rpm		Torque: 1.6 Ncm 300 Rpm		Torque: 1.6 Ncm 300 Rpm		Torque: 1.6 Ncm 300 Rpm		Torque: 1.6 Ncm 300 Rpm		Torque: 1.6 Ncm 300 Rpm		Torque: 1.6 Ncm 300 Rpm		Torque: 1.6 Ncm 300 Rpm		Torque: 1.6 Ncm 300 Rpm		Torque: 1.6 Ncm 300 Rpm		Torque: 1.6 Ncm 300 Rpm									

↑ Excellent Cutting Efficiency

↪ Maximum File Flexibility

EndoArt Smart rotary files have fixed %4 taper allow you to work more safely in narrow canals. The developed special heat treatment technology ensures that the Ni-Ti wire has a much more flexible and durable structure and also turns it into a Gold or Blue color. In addition to its high sharpness efficiency, EndoArt Smart Gold has a higher fracture resistance than standard Ni-Ti files. EndoArt Smart Blue allows pre-bending before entering the canal thanks to its high flexibility efficiency with special control memory and it has a much higher fracture resistance than standard Ni-Ti files. It adapts to the shape of the canal in cases with high slopes.

Smart Kit 15/04-20/04-25/04-30/04-35/04-40/04 19/08-10/04-15/04-20/04-25/04-30/04 19/08-15/04-20/04-25/04-30/04-35/04 19/08-16/02-20/04-25/04-30/04-35/04

**Smart Gold Rotary Files**

- Rotation Files
- Regular Taper
- Better Cutting Efficiency
- x2 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box



**Smart Blue Rotary Files**

- Rotation Files
- Regular Taper
- Better Flexibility
- X4 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box

**Wire** : Gold  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files

**Wire** : Blue  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files







TOUCH

10/06	15/06	20/06	25/06	30/06	35/06	40/06	50/06
Shaping		Finishing					
Torque: 3.0 Ncm 300 Rpm		Torque: 2.6 Ncm 300 Rpm		Torque: 2.6 Ncm 300 Rpm		Torque: 2.6 Ncm 300 Rpm	
Torque: 2.6 Ncm 300 Rpm		Torque: 2.6 Ncm 300 Rpm		Torque: 2.6 Ncm 300 Rpm		Torque: 2.6 Ncm 300 Rpm	

Excellent Cutting Efficiency  
Maximum File Flexibility

10/06	15/06	20/06	25/06	30/06	35/06	40/06	50/06
Shaping		Finishing					
Torque: 3.0 Ncm 300 Rpm		Torque: 1.8 Ncm 300 Rpm		Torque: 1.8 Ncm 300 Rpm		Torque: 1.8 Ncm 300 Rpm	
Torque: 1.8 Ncm 300 Rpm		Torque: 1.8 Ncm 300 Rpm		Torque: 1.8 Ncm 300 Rpm		Torque: 1.8 Ncm 300 Rpm	

EndoArt Touch rotary files fixed %6 taper allow you to prepare at a standard angle along the canal. The developed special heat treatment technology ensures that the Ni-Ti wire has a much more flexible and durable structure and also turns it into a Gold or Blue color. In addition to its high sharpness efficiency, EndoArt Touch Gold has a higher fracture resistance than standard Ni-Ti files. EndoArt Touch Blue allows pre-bending before entering the canal thanks to its high flexibility efficiency with special control memory and compared to standard Ni-Ti files, it has a much higher fracture resistance. It adapts to the shape of the canal in cases with high slopes.

Touch Kit 15/06-20/06-25/06-30/06-35/06-40/06 19/08-15/04-20/04-20/06-25/06-30/06 19/08-10/06-15/06-20/06-25/06-30/06 19/08-15/06-20/06-25/06-30/06-35/06

Touch Gold Rotary Files

- Rotation Files
- Regular Taper
- Better Cutting Efficiency
- x2 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box



Touch Blue Rotary Files

- Rotation Files
- Regular Taper
- Better Flexibility
- X4 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box

Wire : Gold  
Length : 21mm, 25mm, 31mm  
Packaging : 6 Files

Wire : Blue  
Length : 21mm, 25mm, 31mm  
Packaging : 6 Files



WhatsApp  
Get information!

# ENDOART<sup>®</sup>

FINAL TOUCH



**Reciprocation**



EXPERT




EndoArt Expert files work with a reciprocating motion that allows you to perform root canal treatment with a single file. It includes alternative tapers, %6 and %4 angles for single-file systems, and allows you to prepare at a standard angle along the canal. EndoArt Expert files, which can be used with endomotors capable of reciprocating movement, work with a rotation angle adjustment of 150° left and 30° right. The developed special heat treatment technology makes the Ni-Ti wire a much more flexible and durable structure. EndoArt Expert Gold files have higher sharpness and higher fracture resistance than standard Ni-Ti files. EndoArt Expert Blue allows pre-bending before entering the canal thanks to its high flexibility efficiency with special control memory and compared to standard Ni-Ti files, it has a much higher fracture resistance.

Expert Kit R25-R25-R25-R25-R40-R50 R25-R25-R40-R40-R50-R50 R25-R25-R40-R50 R25-R40-R50

**Expert Gold Reciprocation Files**

- Reciprocation Files
- Regular Taper
- Better Cutting Efficiency
- x2 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box



**Wire** : Gold  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files

**Expert Blue Reciprocation Files**


- Reciprocation Files
- Regular Taper
- Better Flexibility
- X4 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box

**Wire** : Blue  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files





PEDO



EndoArt Pedo rotary files are produced specially for pediatric cases and has a length of 18 mm. The improved heat treatment technology of EndoArt Pedo ensures that the Ni-Ti wire has a much more flexible and durable structure. EndoArt Pedo files have higher sharpness and fracture resistance than standard Ni-Ti files.

Pedo Kit 15/06-20/04-25/04 15/06-20/04-30/04

**Pedo Gold Rotary Files**

- Rotation Files
- Specially produced for pedodontic cases
- Regular Taper
- Better Cutting Efficiency
- x2 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box

Control Memory  
Heat Treated Wire

Convex Triangle  
Cross - Section



**Pedo Blue Rotary Files**

- Rotation Files
- Specially produced for pedodontic cases
- Regular Taper
- Better Flexibility
- X4 times better Cyclic Fatigue Resistance than conventional Files
- Excellent Strength
- No bounce back
- Respect Original Anatomy of Root Canal
- Great Fracture Resistance
- Lasered Taper and Size information on the shank
- Autoclavable Box

**Wire** : Gold  
**Length** : 18mm  
**Packaging** : 3 Files

**Wire** : Blue  
**Length** : 18mm  
**Packaging** : 3 Files



WhatsApp  
Get information!

## RETREATMENT



EndoArt Retreatment rotary files developed for repeated root canal treatments. Thanks to its high cutting efficiency, it is used in the process of removing the materials in the canals. The system consists of three instruments, D1 - D2 and D3 which are designed for cleaning the canal filling in the (D1) coronal, (D2) middle, and (D3) apical thirds. While D1 has an active tip at the first entry into the root canal, the tips of continuing instruments are not cutting.

Retreatment Kit D1-D1-D2-D2-D3-D3

### Retreatment Files

Rotation Files  
 Specially produced for retreatment cases  
 Designed to remove filling materials like Gutta Percha  
 x2 times better Cyclic Fatigue Resistance than conventional Files  
 Better Cutting Efficiency  
 Great Fracture Resistance  
 Lasered Taper and Size information on the shank  
 Autoclavable Box

**Wire** : Gold  
**Length** : 16mm, 18mm, 22mm  
**Packaging** : 6 Files



## PATH FILE




Path files are used to create a glide (Glide) path. The use of angled files in combination with hand files before use, reduces the risk of step formation and perforation. The size and taper of the file are 16/02. The ends of the file are rounded to prevent the formation of step and zip formation.

### Path Gold Rotary Files

Rotation Files  
 Specially produced for creating glide path  
 Easily reach apex  
 Better Flexibility  
 x2 times better Cyclic Fatigue Resistance than conventional Files  
 Great Fracture Resistance  
 Lasered Taper and Size information on the shank  
 Autoclavable Box

**Wire** : Gold  
**Length** : 21mm, 25mm, 31mm  
**Packaging** : 6 Files



WhatsApp  
 Get information!



## EndoArt K Files

Codes : 06 08 10 15 20 25 30 35 40 45  
50 55 60 70 80 90 100 110 120 130 140

Ass. 15/40 – Ass. 45/80 – Ass. 90/140

Stainless Steel

Length: 21 mm, 25 mm, 31 mm

6 pieces in a box



## EndoArt H Files

Codes : 06 08 10 15 20 25 30 35 40 45  
50 55 60 70 80 90 100 110 120 130 140

Ass. 15/40 – Ass. 45/80 – Ass. 90/140

Stainless Steel

Length: 21 mm, 25 mm, 31 mm

6 pieces in a box



## EndoArt Spreader

Codes: 15 20 25 30 35 40/ Ass. 15/40

Stainless Steel

Length: 21 mm, 25 mm

6 pieces in a box



## EndoArt C Files

Codes: 06 08 10 12.5 15/ Ass. 06-10

Stainless Steel

Length: 25 mm

6 pieces in a blister package

## EndoArt Endomotor and Apex Locator

Cordless endomotor with apex locator function.

**Properties**

- 10 Memory Settings
- 3 Rotation Mode; Clockwise/Anti Clockwise and Reciprocating
- Max 800 Rpm and 4.0 N.Cm
- Micro USB Charging
- 16:1 Contra Angle
- Auto Start in the root canal
- Auto Reverse
- 4 working mode
- Large OLED Display



CE	16:1	OLED screen	AUTO REVERSE	AUTO START-STOP	Reciprocating Motion
----	------	-------------	--------------	-----------------	----------------------

## EndoArt Apex Locator

High Accuracy Apex Locator

**Properties**

- Advanced Multi Frequency Tech
- Adjustable Reference Position
- Auto Calibration



CE	OLED screen	5th Generation	AUTO CALIBRATION
----	-------------	----------------	------------------

## EndoArt Gutta Percha Points

Gutta Percha root canal filling material for easy fillings. It has variable taper and sizes.

**Properties**

- 02/04/06 Taper Options
- Excellent Radiopacity
- Length Marked
- Hand Rolled



## EndoArt Paper Points

Absorbent Paper Point to dry root canals.

**Properties**

- Highly Absorbent
- Sterile Packaging
- Variety Sizes and Tapers
- Length Marked
- Hand Rolled







## EndoArt Fiber Post

Universal resin glass fiber post system for esthetic and metal free restorations. It has different type of diameters which are 1,2 mm, 1,4 mm, 1,6 mm and 1,8 mm.

### Properties

Radiopaque  
 4 drills available inside of the kit  
 Conical tip  
 High resistance  
 Colour coded  
 It is compatible with adhesive resin cements

### Packaging

EndoArt Fiber Post Kit – 5 pcs. 1,2 mm posts, 5 pcs. 1,4 mm posts, 5 pcs. 1,6 mm posts, 5 pcs. 1,6 mm posts, 4 drills assorted  
 Ø 1.2 - Ø 1.4 - Ø 1.6 - Ø 1.8

EndoArt Fiber Post Refill - 5 pcs. 1,2 mm posts

EndoArt Fiber Post Refill - 5 pcs. 1,4 mm posts

EndoArt Fiber Post Refill - 5 pcs. 1,6 mm posts

EndoArt Fiber Post Refill - 5 pcs. 1,8 mm posts



## RubyCal Calcium Hydroxide Paste With Barium Sulphate

Calcium hydroxide premixed paste non-setting for root canal therapy and vital pulp therapy. It comes with application tips.

**Properties**

- High Radiopacity
- Easy removal
- High pH > 12
- Excellent handling characteristics

**Packaging**

- 2 gr syringe
- 2 gr syringe x 2
- 2 gr syringe x 4



## RubyCal Calcium Hydroxide Powder & Liquid

Calcium hydroxide powder and liquid for root canal therapy

**Properties**

- Temporary Root Filling
- High pH > 12
- Excellent handling characteristics
- Direct and indirect pulp capping
- Base under non-adhesive restorations
- Easy to mix

**Packaging**

- 10 gr powder and 10 ml liquid



## RubyCal Calcium Hydroxide Paste With Iodoform

Calcium hydroxide premixed paste non-setting for root canal therapy and vital pulp therapy. It comes with application tips.

**Properties**

- Radiopacity
- Promote apexification
- Disinfects Canal
- Excellent handling characteristics

**Packaging**

- 2 gr syringe
- 2 gr syringe x 2
- 2 gr syringe x 4



## RubyCal Calcium Hydroxide Powder

Calcium hydroxide powder for root canal therapy.

**Properties**

- Temporary Root Filling
- Direct and indirect pulp capping
- Base under non-adhesive restorations
- Antibacterial Properties
- Easy to mix

**Packaging**

- 30 gr powder



## RubyZoe Zinc Oxide Eugenol Cement

Zinc oxide and eugenol based dental cement for temporary filling.

### Properties

Temporary Filling  
 Root canal sealing  
 Root canal filling  
 Radiopaque  
 Base and lining material under fillings  
 Easy to mix

### Packaging

45gr powder and 15 ml liquid



## RubyZoe Eugenol

Eugenol oil mixing with Zinc oxide powder for temporary filling.

### Properties

%100 Eugenol (Pharma Grade)  
 Easy to mix

### Packaging

15 ml liquid



## RubyZoe Zinc Oxide Powder

Zinc oxide powder mixing with eugenol for temporary filling.

### Properties

%100 Zinc Oxide (Pharma Grade)  
 Base and lining material under fillings  
 Temporary Fixing crown and bridges  
 Indirect pulp capping  
 Easy to mix

### Packaging

45gr powder



## RubyEdta 17% EDTA Gel

It is 17% EDTA based root canal preparation gel for lubricating and cleaning of root canal. It comes with application tips.

**Properties**

- Easily moving
- Antibacterial effect
- Reduce file breakage

**Packaging**

- 3 ml syringe
- 3 ml syringe x 2
- 3 ml syringe x 4



## EndoArt Chlorex Gel

EndoArt Chlorex is %2 chlorhexidine gel for preparation of root canal.

**Properties**

- 2% chlorhexidine digluconate
- Easier root canal treatment
- Minimum risk of root canal breakage

**Packaging**

- 2 gr in a syringe



## EndoArt EDTA Solution

EndoArt EDTA is irrigation solution for preparation and cleaning of root canals, remove smear layer. %17 EDTA or %5 EDTA options are available.

**Properties**

- Clean root canals rapidly
- Neutral PH
- Helps to open calcified canals

**Packaging**

- 250 ml liquid



## EndoArt Chlorex

EndoArt Chlorex is %2 chlorhexidine digluconate irrigation liquid for root canal treatment.

**Properties**

- 2% chlorhexidine digluconate
- Improves mechanical preparation and cleaning of root canals

**Packaging**

- 250 ml liquid



## EndoArt Sodium Hypochlorite

It is 5.25% sodium hypochlorite irrigation liquid for root canal treatment.

### Properties

It has surfactant increase penetration to dentinal tubules.  
Removes smear layer  
Effective canal cleaning  
Strong antibacterial effect

### Packaging

250 ml liquid



## EndoArt Citric Acid

It is %40 citric acid irrigation liquid for root canal treatment.

### Properties

Removes smear layer from the dentinal walls  
Effective canal cleaning  
Cleaning calcium hydroxides from the root canals  
Improves the mechanical cleaning of the root canal

### Packaging

250 ml liquid



## Ruby Iodoform Powder

Iodoform is used for root canal treatment. Can be used for filling of root canals after depulping, chronic periodontitis, wound dressing.

### Properties

Pure Iodoform (Pharma Grade)  
High PH 12,5

### Packaging

30 gr powder



## Ruby Gutta Solver

RubyGuttaSolver is used for removing of gutta percha materials during re-treatment of root canals.

### Properties

%100 Orange Oil  
Non toxic, natural  
Easy to wash out

### Packaging

10 ml liquid

**ruby dent**





incidental



## RubyFill Nano Composite

RubyFill Nano is radiopaque nano hybrid composite system which can be the photocurable applied in the anterior and posterior field.

### Properties

- Universal Nano Composite
- It has light refraction property like a normal tooth due to the nano fillers in it
- It has high filler rate
- It can be polished high with high degree
- It has high resistance to abrasion
- It provides ease of use

### Packaging

- Refill – 4 gr in a syringe
- Kit 1 - 4 gr syringe x 6 pcs. , 5 ml RubyBond x 1 pcs. 3 ml RubyEtch x 1 pcs. and Accessories
- Shades – A1 – A2 – A3 – A3,5 – B1 – B2 – C1 – D3 – OA2 – Extra Bleach



## RubyComp Nano Composite

RubyComp Nano is radiopaque nano hybrid composite system which can be the photocurable applied in the anterior and posterior field.

### Properties

- Higher Esthetic
- Low shrinkage
- It has high filler rate
- It can be polished with high degree
- It has high resistance to abrasion
- It provides ease of use
- Long durability
- Non sticky

### Packaging

- Refill – 4 gr in a syringe
- Kit 1 - 4 gr syringe x 6 pcs. , 5 ml RubyBond x 1 pcs. 3 ml RubyEtch x 1 pcs. and Accessories
- Shades – A1 – A2 – A3 – A3,5 – B1 – B2 – C1 – D3 – OA2 – Extra Bleach



ISO 4049

CLASS 2A



WhatsApp  
Get information!





## RubyComp P Posterior Composite

RubyComp P is micro hybrid composite system which can be the photocurable applied in the posterior field.

### Properties

- Non stick to instruments
- It has surface hardness is like real enamel
- It has high intensity
- It has a perfect consistence can be shaped/formed easily
- It has high filler rate
- It resists to chewing force
- There is not discoloration after the polimeration.

### Packaging

Refill – 4 gr in a syringe  
 Kit 1 - 4 gr syringe x 6 pcs. , 5 ml RubyBond x 1 pcs. 3 ml RubyEtch x 1 pcs. and Accessories  
 Shades – A1 – A2 – A3 – A3,5



## RubyFlow Flowable Composite

RubyFlow is flowable light curing radiopaque composite.

### Properties

- Low shrinkage, low wear
- High gloss polishability
- Beautiful esthetic restorations
- Very low surface roughness
- Excellent handling
- Non sticky
- Very good marginal adaptation

### Packaging

Refill – 2 gr in a syringe  
 Kit 1 - 2 gr in a syringe x 2 pcs.  
 Kit 2 - 2 gr syringe x 4 pcs.  
 Kit 3 - 2 gr syringe x 4 pcs. , 5 ml RubyBond x 1 pcs. 3 ml RubyEtch x 1 pcs. and Accessories  
 Shades – A1 – A2 – A3 – A3,5 – B1 - OA2



ISO 4049

CLASS 2A



WhatsApp  
Get information!



## RubyBond Adhesive System

RubyBond is light curing single component 5. generation bonding agent for the adhesive restorations. It universal and compatible with all composites.

### Properties

Low abrasion  
Elastic connection  
5. generation

### Packaging

5 ml in a bottle  
5 ml in a bottle + 3ml Etching



## RubySEBond Self Etching Adhesive System

RubySEBond is self-etching light curing one component 7 generation adhesive for the adhesive restorations.

### Properties

Very strong bonding to enamel and dentine  
No extra etching with phosphoric acid  
Long lasting restorations  
Time savings  
7. generation

### Packaging

5 ml in a bottle



CLASS 2A



 **WhatsApp**  
Get information!



## RubyEtch Etching Gel

RubyEtch is an etching gel containing 37% phosphoric acid. Delivered in 3 ml. and 50 ml syringe with very thin application needles.

### Properties

Easily Rinse  
Blue coloured  
Perfect Surface Retention

### Packaging

3 ml syringe  
3 ml syringe x 2  
3 ml syringe x 4  
50 ml syringe



## RubyEtch HF Porcelain Etch

RubyEtch HF is 9.5% Buffered Hydrofluoric acid in gel form. Using to improve the mechanical bonding between porcelain and adhesive agent.

### Properties

Porcelain Restoration Repair  
Buffered Formula  
Easily Rinse

### Packaging

3 ml syringe  
Porcelain Repair Kit - 3 ml syringe + Silane



CLASS 2A



WhatsApp  
Get information!



## RubyLiner with Glass Ionomer

Light curing glass ionomer based cavity liner to protect pulp of the tooth. It is used under direct and indirect restorations, like composites, amalgams, ceramics and another metal restorations.

### Properties

Excellent biocompatibility  
Releases flouride  
Radiopaque  
Excellent handling characteristics

### Packaging

2 gr syringe  
2 gr syringe x 2  
2 gr syringe x 4



## RubyCal LC Calcium Hydroxide

Light curing calcium hydroxide for pulp capping and liner material under dental filling.

### Properties

Direct & indirect pulp capping  
Light curing  
Single component

### Packaging

2 gr syringe  
2 gr syringe x 2  
2 gr syringe x 4



ISO 4049

CLASS 2A



WhatsApp  
Get information!



## Ruby Pit and Fissure Sealant

Light curing sealant with fluoride releasing for occlusal surfaces. It can be used as preventive products.

### Properties

Fluoride releasing  
Radiopaque  
Easy application  
Low shrinkage

### Packaging

2 gr syringe  
2 gr syringe x 2  
2 gr syringe x 4



## FlexDam Gingiva Barrier

Light curing resin-based gingiva barrier to protect soft tissue.

### Properties

Easily removable after treatment  
It minimizes heat and soft tissue sensitivity during curing  
Blue shade  
5 second curing

### Packaging

2 gr syringe  
2 gr syringe x 2  
2 gr syringe x 4



CLASS 2A



WhatsApp

Get information!



## Rubydent Curing Light

Rubydent Curing Light offers effective and quality curing performance with ergonomic design and metal head.

- 2 mm resin curing in 1 second
- Long life battery / Over 1300 circles under 10 sec mode
- Stainless steel head
- Replaceable and 360 degree rotating head
- Ergonomic design
- Black or white colour options

Model	: MaxCure7 Power
Input	: 100-240VAC 50Hz/60Hz Power
Output	: 5V DC/1A Wave
length	: 385nm-515nm
Light intensity	: 1000mW/cm2- 2500mW/cm2
Battery	: 2600mA

CE Class I



WhatsApp  
Get information!




## RubyShine Prophylaxis Paste

RubyShine Prophylaxis Paste is specially formulated of polishing and cleaning agents which has nano particuls in it. Available in mint, strawberry, bubblegum and green apple flavor.

### Packaging

30 gr in a bottle  
150 gr in a bottle  
340 gr in a bottle



Class I



WhatsApp  
Get information!



## RubyCap Capsule Amalgam

High copper alloy which contains a mixture of spherical and irregular particles. Predosed in self activating capsules. Gamma-2 free. It has 4 spill options as spill 1, spill 2, spill 3 and spill 5.

### Indications

Stress bearing fillings in posterior teeth (Class I and II)  
Retrograde root canal filling after apex resection

### Features

Mixture of spherical and irregular particles  
High copper alloy  
High affinity with mercury  
Firm condensation  
Smooth carving  
Excellent condensability

### Composition

45% Ag – 30,5% Sn – 24% Cu – 0,5% Zn  
45% Ag – 31% Sn – 24% Cu (Zinc Free)

### Packaging

50 capsule in a jar



## RubyCap NG Capsule Amalgam

High copper alloy which contains a mixture of spherical and irregular particles. It contains higher silver (55% Ag) than RubyCap. Predosed in self activating capsules. Gamma-2 free. It has 4 spill options as spill 1, spill 2, spill 3 and spill 5.

### Indications

Stress bearing fillings in posterior teeth (Class I and II)  
Retrograde root canal filling after apex resection

### Features

High silver / high copper alloy  
Easy to use

### Composition

55% Ag – 29% Sn – 15,8% Cu – 0,2% Zn

### Packaging

50 capsule in a jar



ISO 20749

CLASS 2A



WhatsApp  
Get information!





## RubyCap High Capsule Amalgam

High silver alloy which contains a mixture of spherical and irregular particles. Predosed in self activating capsules. Gamma-2 free. It has 4 spill options as spill 1, spill 2, spill 3 and spill 5.

### Indications

Stress bearing fillings in posterior teeth (Class I and II)  
Retrograde root canal filling after apex resection

### Features

Mixture of spherical and irregular particles  
High silver content  
Exceptionally durable  
Easy to process  
Relatively fast hardening

### Composition

%69,2 Ag %18,6 Sn %11,9 Cu %0,3 Zn  
%69,2 Ag %18,9 Sn %11,9 Cu (Zinc Free)

### Packaging

50 capsule in a jar

### Technical Properties

A well balanced alloy for quality fillings;  
Spill 1 – 400 mg Alloy / 360 mg Mercury  
Spill 2 – 600 mg Alloy / 540 mg Mercury  
Spill 3 – 800 mg Alloy / 720 mg Mercury  
Spill 5 – 1200 mg Alloy / 1080 mg Mercury

Working Time: 4 minutes  
Carving & burnishing: 5 – 9 minutes

Conforms to:  
ISO 20749 Pre-capsulated Dental Amalgam

Compressive Strength;  
1h < 397 MPa  
24h < 437 MPa

Creep : %0,99  
Dimensional Change During Hardening : %0,01

RubyCap I RubyCap II RubyCap III



RubyCap NG I RubyCap NG II RubyCap NG III



RubyCap High I RubyCap High II RubyCap High III



ISO 20749

CLASS EA





## RubyGam Powder Amalgam

High copper alloy which contains a mixture of spherical and irregular particles. Gamma-2 free.

### Indications

Stress-bearing fillings in posterior teeth (Class I and II)  
Retrograde root canal filling after apex resection

### Features

Mixture of spherical and irregular particles  
High copper alloy  
High affinity with mercury  
Firm condensation  
Smooth carving  
Excellent condensability

### Composition

45% Ag – 30,5% Sn – 24% Cu – 0,5% Zn

### Packaging

Boxes of 30 gr

## RubyGam High Powder Amalgam

High silver alloy which contains a mixture of spherical and irregular particles. Gamma-2 free.

### Indications

Stress bearing fillings in posterior teeth (Class I and II)  
Retrograde root canal filling after apex resection

### Features

A mixture of spherical and irregular particles  
High silver content  
Exceptionally durable  
Easy to process  
Relatively fast hardening

### Composition

69,2% Ag – 18,6% Sn – 11,9% Cu – 0,3% Zn

### Packaging

Boxes of 30 gr



**ruby | dent**



## Ruby Cap

STABILITY

HIGH STRENGTH

EXCELLENT HANDLING

NON GAMMA 2

## RubyPlaton Polishing Discs

Has different kind of discs which has different grins and different sizes. It allows excellent results in polishing. Flexible and easy to use


### Properties

14 mm blue, yellow, white, green, brown, pink

- **White** – Extra Fine for Polishing
- **Yellow** – Fine for Finishing
- **Green** – Medium for Countoring
- **Blue** – Coarse for Gross Reduction
- **Brown** – X Coarse for Gross Reduction
- **Pink** – XX Coarse for Gross Reduction

### Packaging

80 pcs discs and 1 mandren in a plastic box  
240 pcs discs and 4 mandren in a plastic box



## RubyPlaton Polishing Strips

Has different kind of strips which has different grin sizes. They are produced combining two different types of grins.

### Properties

It allows excellent results in polishing. Flexible and easy to use

- **White** – Extra Fine for Polishing
- **Yellow** – Fine for Finishing
- **Green** – Medium for Countoring
- **Blue** – Coarse for Gross Reduction

### Packaging

75 pcs strips in a plastic box  
80 pcs discs and 1 mandren + 75 pcs strips in a plastic box  
240 pcs discs and 4 mandren + 75 pcs strips in a plastic box





## RubyPlaton Wooden Wedges

RubyPlaton wooden wedge it has anatomic form for you to make excellent filling.

### Properties

- Superthin short
- Thin short
- Thin long
- Medium short
- Medium long

### Packaging

100 pcs. wedges in a plastic box  
 200 pcs. wedges in a plastic box  
 500 pcs. wedges in a plastic box



## RubyPlaton Transparent Wedges

RubyPlaton transparent wedge is formed by special transparent polymer for excellent filling.

### Properties

Not easy to break long lived with high quality and durability  
 Easy to use  
 Anatomic form  
 Small size - Medium size

### Packaging

100 pcs in a plastic box

CE Class I



WhatsApp  
 Get information!



## RubyMark Articulating Paper Straight

Great sensitive articulating paper to get simple marking. It is straight form on the 40µ, 80µ and 200µ.

### Properties

Non sticky  
Smooth surface  
Double sided options.

### Packaging

40µ – Blue/Blue Straight 144 pcs. paper  
40µ – Red/Red Straight 144 pcs. paper  
40µ – Blue/Red Straight 144 pcs. paper  
80µ – Blue/Blue Straight 144 pcs. paper  
80µ – Red/Red Straight 144 pcs. paper  
80µ – Blue/Red Straight 72 pcs. paper  
80µ – Blue/Red Straight 144 pcs. paper



Class I



WhatsApp  
Get information!



## RubyMark Articulating Paper Horseshoe

Great sensitive articulating paper to get simple marking. It is horseshoe form on the 80 $\mu$ .

### Properties

Non sticky  
Smooth surface  
Double sided options.

### Packaging

80 $\mu$  – Blue/Red Horseshoe 72 pcs. paper

## RubyMark Articulating Paper 200 Micron Straight

Great sensitive articulating paper with special fleece paper to get simple marking. It is straight form on the 200 $\mu$ .

### Properties

Non sticky  
Smooth surface  
Most Accurate Occlusion

### Packaging

200 $\mu$  – Blue Straight 300 pcs. paper  
200 $\mu$  – Red Straight 300 pcs. paper



CE Class I



WhatsApp  
Get information!



## RubyCrownRemover

It is practical crown remover for easy procedure. It removes crown using the way of impact method with arc.

### Properties

- Can be sterilized at 134 °C
- 3 types of remover tips
- Different types of impact



## RubyPlaton Cotton Roll

High quality cotton rolls made from %100 cotton. They are non-linting, non-sterile, soft texture rolls.

### Properties

- It retracts cheek effective
- It does not adhere to mucous membrane
- Highly absorbent
- Size: 4 x 1 cm
- Maintain shape when wet

### Packaging

- No: 1 / 500 pieces bag
- No: 2 / 500 pieces bag
- No: 3 / 500 pieces bag

CE Class I



WhatsApp  
Get information!





## Arma Dental Study Model


Dental study model for students. It shows great and real-like experience during application.

### Properties

Hard Plastic jaw base  
 Non-breakable  
 The students thus practice under "human" conditions  
 Custom made screws for tooth fixing  
 Natural anatomy and occlusion

### Packaging

1 piece dental study model with screwed teeth



## Arma Dental Phantom Teeth

Dental study model teeth for students. It shows great and real-like experience during application. Teeth are made from melamine material.

### Properties

Custom made screws for tooth fixing  
 Natural anatomy and occlusion

### Packaging

1 set dental study model teeth and screws with teeth chart

# ENDOART<sup>®</sup>

FINAL TOUCH



 incidental



For more information;  
export@incidental.com.tr



**+905302252254**



**incidental**

**İnci Dental Tıbbi Malz. San. ve Tic. Ltd. Şti.**

Defterdar Mah. Otakçılar Cad.

E Blok Sit. (Has On Ticaret Merkezi) No:74/4 Eyüpsultan / İstanbul / Turkey

Tel : +90 212 632 63 45 Fax : +90 212 632 97 56

export@incidental.com.tr www.incidental.com.tr